

# SteadyDraw Custom Size Request

Want a SteadyDraw for a syringe or vial size we don't make yet? Fill this out and I'll build it.

## 1 — WHAT ARE YOU USING?

Syringe brand & model (e.g., BD Ultra-Fine, EasyTouch): \_\_\_\_\_

Syringe capacity (mL / cc): \_\_\_\_\_

Needle: fixed (insulin-style) or detachable (Luer)? \_\_\_\_\_

Vial capacity (mL): \_\_\_\_\_

Vial brand / medication packaging (if known): \_\_\_\_\_

## 2 — SYRINGE MEASUREMENTS (diagram on page 2)

- A** Barrel outer diameter \_\_\_\_\_ mm
- B** Barrel length — needle hub to finger flange \_\_\_\_\_ mm
- C** Finger flange width — wing tip to wing tip \_\_\_\_\_ mm
- D** Plunger thumb-rest width \_\_\_\_\_ mm
- E** Needle length (if attached) \_\_\_\_\_ mm
- F** Shoulder width — where the barrel steps down \_\_\_\_\_ mm

## 3 — VIAL MEASUREMENTS (diagram on page 3)

- A** Total height, cap included \_\_\_\_\_ mm
- B** Body diameter (widest part) \_\_\_\_\_ mm
- C** Cap diameter \_\_\_\_\_ mm


## HOW TO MEASURE && SEND

- Measure in millimeters (mm) — calipers are best, a ruler works fine.
- Syringe: plunger pushed all the way in. Vial: cap on.
- Photos of the syringe and vial next to a ruler help enormously — send those even if you can't measure.
- Reply with the letters and your numbers (example: Syringe A = 6.4 mm, B = 54 mm ...)
- Send via Etsy message, Makerworld comment, or email: Samuelcon41@gmail.com

## Page 2 — Syringe Dimensions

### Custom Size Request — SYRINGE

Please measure in millimeters (mm) with the plunger pushed all the way in




- |                                                      |                                                   |
|------------------------------------------------------|---------------------------------------------------|
| <b>A</b> Barrel outer diameter (across the tube)     | <b>D</b> Plunger thumb-rest width                 |
| <b>B</b> Barrel length — needle hub to finger flange | <b>E</b> Needle length (if attached)              |
| <b>C</b> Finger flange width — wing tip to wing tip  | <b>F</b> Shoulder width — where barrel steps down |
- + Capacity (ml/cc) & brand if known

A photo of the syringe next to a ruler helps too! — AL & Oats 3D

## Custom Size Request — VIAL

Please measure in millimeters (mm), cap on



A photo of the vial next to a ruler helps too! — AL & Oats 3D